

# Bordered Foam

Adhesive Bordered Hydrophilic Foam Dressing – Waterproof



**Uses:**

- Manages low to moderate exuding wounds
- Chronic & acute wounds
- Partial & full thickness wounds
- Diabetic ulcers
- Donor sites
- Pressure ulcers – Stages II, III, IV

**Comparable to \*OptiFoam Adhesive Island Dressing**

- ▶ Fluid handling system
- ▶ Exudate absorbed by dressing
- ▶ Waterproof outer layer
- ▶ Conforms to difficult wound site
- ▶ Contraindicated for third degree burns
- ▶ Minimizes the risk of per wound skin maceration and erosion

**LATEX FREE**

Catalog #	Description	Per Box	NDC # 61924
00298E	BorderedFoam 4" x 4" Dressing	10	298-44

\*Medline

**DermaRite**  
INTELLIGENT SKIN & WOUND CARE  
DermaRite Industries LLC  
800-337-6296 • www.dermarite.com  
800-DERMA96 • info@dermarite.com