

# Bordered Gauze <sup>TM</sup>

## Absorbent Island Dressing

Dressing changes will depend on condition of wound and the amount of exudate



### Uses:

- To manage surgical incisions
- Superficial lacerations, abrasions
- Stage I-IV pressure ulcers
- Partial and full thickness wounds
- Facilitates wicking of fluid
- 1<sup>st</sup> and 2<sup>nd</sup> degree burns

### Comparable to \*Steri-Pad

- ▶ Protects the wound surface
- ▶ Soft and flexible
- ▶ Helps to conform dressing to difficult wound site
- ▶ Gentle on skin
- ▶ Non-Adherent contact layer permeable
- ▶ An absorptive pad
- ▶ Adhesive bordered non-woven cloth tape
- ▶ Absorptive dressing consists of three layers

### LATEX FREE

Catalog #	Description	Per Case	NDC # 61924-
00255	Bordered Gauze 4x4 (2x2 pad)	100 per bag	255-44
00256	Bordered Gauze 6x6 (4x4 pad)	100 per bag	256-66



DermaRite Industries LLC

800-337-6296 • [www.dermarite.com](http://www.dermarite.com)

800-DERMA96 • [info@dermarite.com](mailto:info@dermarite.com)

\*Johnson & Johnson