

DermaGinate™

Calcium Alginate Wound Dressing



Uses:

- A Topical Wound Dressing For Moderate To High Exudating Wounds

- ▶ Conforms readily to wound bed
- ▶ Will maintain hydrophilic environment
- ▶ Creates a moist wound environment
- ▶ When gelled, DermaGinate assists in maintaining moist wound healing environment
- ▶ May absorb up to 17 times its own weight
- ▶ Derived from natural seaweed
- ▶ Will not stick to healing tissue

LATEX FREE

Catalog #	Description	Per Box	NDC # 61924-
00270E	DermaGinate 4x4	10	270-44
00275E	DermaGinate 12" Rope	5	275-12



DermaRite Industries LLC
 800-337-6296 • www.dermarite.com
 800-DERMA96 • info@dermarite.com