

HydraFoam™

Hydrophilic Foam Dressing



Uses:

- Manage chronic & acute wounds
- Lower extremity ulcers, pressure ulcers (stages 2,3,4), donor sites, partial and full thickness wounds
- Low to moderate exuding wounds
- Partial & full thickness wounds

Comparable to *OptiFoam **Hydrasorb

- ▶ Maintains a moist wound environment.
- ▶ Prevents bacterial contamination.
- ▶ Protects and cushions wound.
- ▶ Conforms to awkward to dress areas.
- ▶ Absorbs up to 20 times its own weight.

LATEX FREE

Catalog #	Description	Per Box	NDC# 61924
00294E	Hydrafoam 2x2 Dressing	10	294-22
00295E	Hydrafoam 4x4.25 Dressing	10	295-44
00296E	Hydrafoam 6x6 Dressing	10	296-66

*Medline

**Kendall Healthcare Products



DermaRite Industries LLC

800-337-6296 • www.dermarite.com

800-DERMA96 • info@dermarite.com