


SilverDerm7

Antimicrobial Dressing / Ionic Silver Release

**Advanced
7 Day
Silver
Dressing**



- Uses:**
- Chronic wounds
 - Burns
 - Traumatic wounds
 - Surgical wounds
 - Leg ulcers
 - Diabetic foot ulcers

Equivalent to *Silverlon

- ▶ Local broad spectrum antimicrobial activity
- ▶ Sustained delivery of Silver ions
- ▶ Within 30 minutes Silver ions released
- ▶ Hypo-Allergenic
- ▶ Pain reduction
- ▶ Easy Application
- ▶ Assist with autolysis of necrotic wounds

LATEX FREE

Catalog #	Description	Per Case	NDC # 61924-
00288	Silverderm 7 2x2	10	288-22
00289	Silverderm 7 4x4	10	289-44



*Argentum Medical

DermaRite Industries LLC
800-337-6296 • www.dermarite.com
800-DERMA96 • info@dermarite.com

Zones of Inhibition

Surgical Site Pathogen	Percentage of Isolates 1990-1996 (N=17,671)*	Silverderm 7 Zone of Inhibition**		Silverderm 7 Bactericidal activity
		Number Isolates	Zone Size (mm)	
<i>S. aureus</i>	20%	48	12.5+/-1.1	Yes
<i>Coagulase-Neg Staphylococci</i>	14%	26	12.7+/-1.3	Yes
<i>Enterococcus spp</i>	12%	32	11.7+/-1.5	Yes
<i>E. coli</i>	8%	21	12.7+/-2.6	Yes
<i>P. aeruginosa</i>	8%	33	15.3+/-2.3	Yes
<i>Enterobacter spp</i>	7%	32	11.7+/-1.5	Yes
<i>P. mirabilis</i>	3%	10	10.4+/-1.3	Yes
<i>K. pneumoniae</i>	3%	11	10.7+/-0.9	Yes
<i>C. albicans</i>	3%	1	12.5+/-1.0	Yes
<i>Streptococci spp</i>	4%	34	12.5+/-1.0	Yes

*Data from National Nosocomial Infections Surveillance, CDC, US DHHS

** Data from Heggors JP, Journal of Burns